

FLOOR SENSOR

Floor sensor may be placed at the side of a bed or in a doorway. Floor sensor alerts caregivers when a resident attempts to depart. As soon as floor sensor mat is stepped on, it alerts to caregivers to assist the resident immediately before he/she leaves the room.



correspond to:
NC-2000
NC-100
Gurad System II





Nurse call connection set

Floor sensor Nurse call NC-100 set

NC-100 allows to establish wireless connection between sensor and nurse call system.




* Please refer to page 23 for NC-100



Floor sensor	Nurse call connection set		Nurse call NC-100 set	
	Item Number	Set item	Item Number	Set item
M size set 	EFMS-14/M/Color	Floor mat M size x 1	NC-100F EFMS-14/M/Color	Floor mat M size x 1 NC-100 x 1 set
L size set 	EFMS-14/L/Color	Floor mat L size x 1	NC-100F EFMS-14/L/Color	Floor mat L size x 1 NC-100 x 1 set

Floor sensor My MelodyS set


My MelodyS alerts caregivers when a resident steps on the floor mat. My MelodyS has 8 buzzer/melody mode setting. Most suitable for institutions without nurse call facilities.



* Please refer to page 21 for My MelodyS.

Floor sensor	My MelodyS set	
	Item Number	Set item
M size set 	EXC-02S-04/M/color	Floor mat M size x 1 MymelodyS x 1 set
L size set 	EXC-02S-04/M/color	Floor mat L size x 1 MymelodyS x 1 set

* Please choose color between Black and Grey



Floor sensor mat is water resistance, durable and has wide sensor area.

